

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/225186282>

Automated scoring of AFLPs using rawgeno v 2.0, a free R CRAN library

Article in *Methods in molecular biology* (Clifton, N.J.) · June 2012

DOI: 10.1007/978-1-61779-870-2_10 · Source: PubMed

CITATIONS

31

READS

509

3 authors, including:



Nils Arrigo

University of Lausanne

179 PUBLICATIONS 1,806 CITATIONS

[SEE PROFILE](#)



Nadir Alvarez

Natural History Museum of Geneva

384 PUBLICATIONS 3,403 CITATIONS

[SEE PROFILE](#)

Some of the authors of this publication are also working on these related projects:



ECOGEN - Ecosystem change and species persistence over time: a genome-based approach [View project](#)



Selaginella [View project](#)

Automated scoring of AFLPs using RawGeno v 2.0, a free R CRAN

library

Authors: Arrigo Nils^{1,2}, Holderegger Rolf³, Alvarez Nadir²

¹Laboratory of Evolutionary Botany, Institute of Biology, University of Neuchâtel, 11 rue

Emile-Argand, 2009 Neuchâtel, Switzerland

²Department of Ecology and Evolution, University of Lausanne – UNIL Sorge, Biophore

Building, 1015 Lausanne, Switzerland

³WSL Swiss Federal Research Institute, Zürcherstrasse 111, 8903 Birmensdorf, Switzerland

Keywords: Scoring optimization, command lines, graphical user interface, data mining, bin editing

Abstract

Amplified Fragment Length Polymorphisms (AFLPs) are a cheap and efficient protocol for generating large sets of genetic markers. This technique has become increasingly used during the last decade in various fields of biology, including for instance population genomics, phylogeography and genome mapping. Here, we present RawGeno, an R library dedicated to the automated scoring of AFLPs (i.e. the coding of electropherogram signals into ready-to-use datasets). Our program includes a complete suite of tools for binning, editing, vizualizing and exporting results obtained from AFLP experiments. RawGeno can either be used with command lines and program analysis routines or through a user-friendly graphical user interface. We describe

27 the whole RawGeno pipeline along with recommendations for (i) setting the analysis of
28 electropherograms in combination with PeakScanner, a program freely distributed by Applied
29 Biosystems, (ii) performing quality checks, (iii) defining bins and proceeding to scoring, (iv)
30 filtering non-optimal bins and (v) exporting results in different formats.

31

32 **1. Introduction**

33 The amplified fragment length polymorphism (AFLP) technique is increasingly used in
34 phylogeographic and population genomics studies, particularly in non-model organisms for which
35 no prior DNA sequence information is available **(1)**. This relatively cheap technique is based on
36 complete endonuclease restriction digestion of total genomic DNA followed by selective PCR
37 amplification and electrophoresis of a subset of fragments, resulting in a unique, (theoretically)
38 reproducible fingerprint for each individual. Although the AFLP technique is able to generate a
39 large number of informative markers, the success of this method is compromised by different
40 factors **(2)**. For instance, manual scoring relies on visual inspection and subjective interpretation of
41 the electropherotic profiles during a time-consuming and tedious task. In the last decades, several
42 improvements in automatic scoring have been proposed and implemented in commercial software
43 [see **(3)** for a review]. However, until recently, no free open-source software was available to
44 process AFLP data from raw data to ready-to-use presence/absence binary matrices. Two years ago,
45 we developed RawGeno 1.0 **(4)**, a library performing automated binning, scoring and data mining
46 analyses under the R CRAN environment, based on outputs from the freely available
47 electropherogram-analyzing software PeakScanner (Applied Biosystems, Foster City, USA,
48 <http://www.appliedbiosystems.com/peakscanner>). Implementing RawGeno 1.0 solutions in a free
49 environment has provided an accessible and accurate solution to many users (with 415 downloads
50 from <http://sourceforge.net/projects/rawgeno> one year after release). Here, we present RawGeno
51 2.0, an updated version of this software, built around an optimized, less time-consuming algorithm,
52 implementing new features such as binning editing. We provide examples both in a user-friendly
53 and in a command-line interface. In order to allow users to customize queries (which may vary

54 depending on dataset quality/size), we further provide tips for setting analyses and improving output
55 robustness. All stages of a whole analysis are detailed according to the following five sections: (i)
56 importing features; (ii) quality check; (iii) binning and scoring algorithms; (iv) bin filtering; (v)
57 exporting options.

58

59 **2. Program Utilization**

60 ***2.1. Overview of the analysis***

61 Analysing AFLP electropherograms is achieved using two programs: PeakScanner and
62 RawGeno. Whereas PeakScanner detects AFLP peaks along electropherograms and calculates their
63 intensity and size in base pairs (by relying on an internal size standard included in electrophoresis),
64 RawGeno proceeds to the binning and scoring of AFLP electropherograms.

65 RawGeno includes several filters to assess the quality of electropherograms and checks the
66 consistency of binning. In addition, several preliminary analyses are available to the user for
67 making biological inferences and/or remove outlier samples. Finally, several functions allow
68 exporting resulting data sets into properly formatted files for further analyses.

69

70 ***2.2. Building up a scoring project***

71 In RawGeno 2.0 the AFLP scoring project should be organized according to the following
72 procedure.

- 73 a. Create a folder (hereafter “project folder”) from which the project will be managed.
- 74 b. In this folder, add a sub-directory including all electropherogram files (*.fsa). Also add an R
75 shortcut (Windows users) to conveniently launch RawGeno scoring sessions. Right-clicking
76 on this R shortcut allows defining the default working directory of R by specifying it into the
77 “Start into” addressing field of the shortcut. Copy-pasting the project folder address into this
78 field will set-up the working directory of R accordingly.
- 79 c. Create a text-tabulated table listing individuals included in the project (hereafter referred to
80 as “info table”). The info table is optional as the minimal RawGeno analysis can proceed

without it. However, RawGeno includes several functions relying on this table, for instance to label individuals during preliminary analyses or facilitate the sorting and selection of individuals (for example, according to populations or species) during the production of exports. Therefore, the info table should include any additional relevant information that the user would like to consider. It must contain at least the name of individuals (i.e. in a column named “Tag”) and any supplementary information in extra columns.

2.3. Obtaining the raw data from *.fsa files using PeakScanner

The analysis of AFLPs starts by using PeakScanner in order to detect peaks along electropherograms and to calculate their size. The procedure is highly automated, letting the user set peak detection parameters and check the quality of electropherograms. The peak detection parameters are set up using a so-called “Analysis Method”, which is available from the graphical interface in PeakScanner (menu “Resources/Manage Analysis Methods”). Typically, a proper peak detection attempts to detect only peaks that are biologically relevant and exclude peaks only reflecting technical background noise.

We advise to set the “Analysis Method” using the following guidelines.

- a. Prior to the detection of peaks *per se*, a light smoothing of electropherograms might be desirable, in order to eliminate small secondary peaks due to technical background noise.
- b. The detection of peaks is achieved through a “sliding window” analysis that inspects electropherograms locally. Within the inspected region, PeakScanner first creates a modelled version of the electropherogram by fitting a polynomial curve to the data. Peaks are detected according to this modelled signal, based on their absolute width. Therefore, the detection sensitivity is adjusted by modifying the width of the sliding window (i.e. in terms of data points, the smaller it is, the more sensitive the procedure becomes), the goodness of fit reachable by the polynomial curve (again, increasing the polynomial degree of fitting increases the detection sensitivity) and the minimal width above which a peak is recorded as present. We advise to use the default parameters as a starting point, as they have been shown

to provide reliable results (4, 5): set 15 points for the sliding window width, use a third degree polynomial curve and consider peaks that at least have two measurement points of half-width.

- c. Downstream to peak detection, PeakScanner filters peaks according to their absolute fluorescence intensity, i.e. the peak height, measured in relative fluorescent units (rfu). Visually checking electropherograms obtained from blank samples generally helps to adjust the fluorescence threshold to the upper limit of the technical background noise. While some applications might benefit from considering only peaks with a strong fluorescence (e.g. greater than 150 rfu to provide conservative estimates for band presence statistics), most users will prefer using a more permissive threshold at this stage and apply *a posteriori* filtering strategies based on bin quality statistics (6, 7). We advise to use 50 rfu as a minimal fluorescence for considering individual AFLP peaks.
- d. Save the customized “Analysis Method” in order to use it during electropherogram analysis.
- e. Once the “Analysis Method” is set up, import the electropherograms (stored as *.fsa files) into PeakScanner using the “Add Files” button.
- f. Define the size standard and the “Analysis Method” to be used for all individuals included in the project (set this information for the first individual, then select the columns “Size Standards” and “Analysis Method” and use the “ctrl+D” keyboard shortcut to apply these settings to the remaining individuals).

The detection and sizing of peaks is processed using the “Analysis” menu, from the graphical interface. Once achieved, electropherograms can be visualized and compared among individuals. This might help identifying AFLP reactions that were not successful (e.g. individuals with a systematically low fluorescence or showing abnormal peaks). Removing such individuals prior to the RawGeno analysis will help to enhance the final quality of scoring.

The PeakScanner analysis ends with a simple export process in which the list of peaks detected throughout the complete set of analyzed individuals is stored in a table. This is achieved using the “Export/Export Combined Table” menu, producing a text-tabulated file containing the size, height,

135 area and width of all detected peaks (this can be checked using the “Edit Table Settings” menu).

136

137 ***2.4. From raw data to ready-to-use matrices using RawGeno***

138 *2.4.1. Installing RawGeno*

139 RawGeno is freely available from <http://sourceforge.net/projects/rawgeno> as a zip file. In
140 windows, the installation is achieved either using the graphical user interface of R (menu
141 “Packages/Install package(s) from local zip files”) or the following command line in the R console:

```
142 utils:::menuInstallLocal()
```

143

144 Installing RawGeno with Linux requires decompressing RawGeno.zip into the library folder of
145 R. Using the shell command line, this is done as follows:

```
146 sudo mv RawGeno.zip /usr/lib64/R/library/RawGeno.zip
```

```
147 cd /usr/lib64/R/library/RawGeno.zip
```

```
148 unzip RawGeno.zip
```

```
149 rm RawGeno.zip
```

150

151 Finally, RawGeno requires the installation of two companion packages: *vegan* and *tkrplot*, that
152 both are available from usual R CRAN repositories (Linux users should see Note 1). Their
153 installation is achieved either using the graphical user interface of R (menu “Packages/Install
154 package(s) from CRAN”) or with the following command lines (prompted into the R console):

```
155 install.packages("vegan")
```

```
156 install.packages("tkrplot")
```

157

158 *2.4.2. Importing PeakScanner results*

159 All following steps are performed in the R CRAN environment, using the R shortcut described
160 above (or ensuring that the correct working directory has been selected). Once the RawGeno
161 package has been installed, it can be called applying the following command line:

162 a. Call RawGeno, vegan and tkrplot as libraries into R and launch the graphical user interface
163 `require(RawGeno)`
164 `require(vegan)`
165 `require(tkrplot)`
166 `RawGeno()`

167

168 b. Importing the PeakScanner text-tabulated file in RawGeno can then be done using the
169 graphical user interface (menu “RawGeno/1. Files/Electroph./PeakScanner (*.txt)”) or using
170 the following command lines:

171 Choose interactively the PeakScanner file

172 `myfile=tk_choose.files(caption='Choose PeakScanner File')`

173 `OPENAFLP(myfile, pksc=T, dye="B")`

174

175 Or explicitly specify the path of file of interest

176 `mypath="C:/MyDocuments/MyPeakScannerFile.txt"`

177 `OPENAFLP(mypath, pksc=T, dye="B")`

178

179 During importation, RawGeno handles a single dye color at a time, which is user-specified and
180 considers the “dye” parameter with the following values: "B" (blue; FAM), "G" (green; HEX), "Y"
181 (yellow; NED), “R” (red; ROX) or “O” (orange; LIZ). If electrophoresis was achieved using several
182 dyes simultaneously (e.g. multiplexing of PCR products), each dye must be analysed separately in
183 RawGeno. Datasets obtained from several dyes can be merged *a posteriori* in a final binary table
184 (see below).

185

186 2.4.3. Quality check

187 Because the detection of AFLP peaks is based on a defined threshold, it is not easy to handle
188 reactions showing electropherograms with varying intensities (see Note 2). When improperly

189 handled, such a situation leads to the inclusion of samples characterised by many false-absences in
190 the final dataset. Although the only way to correctly address this issue is a robust wet-lab protocol,
191 RawGeno still attempts limiting the influence of low quality AFLPs on binning and scoring, by
192 filtering individuals that were unsuccessful. Here, the variability in the number of peaks detected
193 per individual is used as a proxy of AFLP reactions quality. Empirical evidence shows that this
194 statistics is dependent of the specific dataset used and the biological organism studied. The lower
195 bound of this distribution most generally includes individuals with low AFLP intensities,
196 characterized by many AFLP peaks that remain undetected in the electropherograms. Because such
197 individuals usually represent a small fraction of the complete project, we advise removing them
198 from the dataset. The upper bound of the distribution can either reflect a biologically relevant signal
199 (e.g. hybridization) or a technical bias (e.g. contamination, odd PCR reaction). Such individuals
200 should be either discarded or identified as outliers for proper interpretation in further analyses.

201 RawGeno includes the two following options for checking quality of individuals:

- 202 a. Filter individuals according to the number of detected peaks, by manually selecting
203 individuals that should be kept for further analysis or by using a dedicated device (Fig 1 a).
204 The command line version of this operation relies on percentiles and conserves individuals
205 within 5%-95% bounds of the detected peaks distribution.

206 Retrieve imported electropherograms

207 `data.electroph=AFLP$all.dat`

208

209 Compute number of AFLP peaks per individual

210 `pk.per.smp=table(data.electroph$sample.id)`

211

212 Compute 5% and 95% quantiles

213 `qtles=quantile(pk.per.smp, probs=c(0.05, 0.95))`

214

215 Determine what samples can be kept

```

216     to.keep=which(pk.per.smp>=qtles[1] & pk.per.smp<=qtles[2])
217
218     Filter electropherograms, keep only retained individuals and update individuals indexing
219     accordingly
220     smp.ok=match(data.electroph$sample.id, to.keep)
221     data.clean=data.electroph[is.na(smp.ok)==F, ]
222     data.clean$sample.id=as.factor(data.clean$sample.id)
223     levels(data.clean$sample.id)=1:length(to.keep)
224     smp.names=AFLP$samples.names
225     smp.clean=smp.names[to.keep]
226     AFLP$all.dat=data.clean
227     AFLP$samples.names=smp.clean
228
229     b. Check the quality of individuals, by taking into account their position in PCR plates (only
230     available from the graphical interface, Fig 1 b). This display helps to highlight systematic
231     biases having a technical origin (e.g. pipetting errors or thermocycler bias) and to identify
232     batches of individuals that were not successful.

```

233

234 2.4.4. Binning

235 Building a presence/absence matrix requires recognizing which AFLP peaks are homologous
 236 across individuals. This procedure relies on the size of peaks along electropherograms and assumes
 237 homology for peaks sharing identical sizes. Because peak sizes are determined empirically using
 238 electrophoresis, measurements generally include technical variations preventing the observation of
 239 strictly identical sizes across homologous AFLP peaks. Indeed, Holland *et al.* (8) reported
 240 measurement variations ranging between 0 and 0.66 bp (with 0.08 bp in average) for replicated
 241 AFLP peaks. Therefore, properly recording the signals of AFLP peaks asks for taking into account
 242 size variations by defining size categories (i.e. “bins”) into which the presence / absence of AFLP

243 peaks are recorded. Bins are characterised by their position along the electropherogram (i.e. the
244 average size of peaks they include) and their width (the size difference between the longest and
245 shortest peaks included in the bin). RawGeno uses a binning algorithm relying on the size of AFLP
246 peaks over all individuals included in the project and defines bins in a way that respects the two
247 following conditions (Fig 2 a).

- 248 a. The first constrain is a maximal bin width. This prevents the definition of too large bins that
249 could lead to homoplasy (i.e. erroneously assigning non-homologous AFLP peaks within the
250 same bin). This limit is set using the MaxBin parameter. We advise to use MaxBin values
251 ranging between 1.5 and 2 bp (see Note 3). Using small values (i.e. MaxBin < 0.5 bp) should
252 be avoided as this generally causes oversplitting, a situation where the presence / absence of
253 homologous AFLP bands are coded using an exaggerated number of bins.
- 254 b. The second constrain prevents the assignment of more than one peak from the same
255 individual within the same bin [i.e. “technical homoplasy”, as defined in (4)]. If such a
256 situation occurs, RawGeno defines two separate bins in which the two peaks are assigned.
257 This constrain can be relaxed by increasing the “MinBin” parameter in order to include two
258 peaks of the same individual in the same bin (when not exceeding the “MinBin” value in size
259 difference). Such a relaxing might be desirable, for instance when artefactual peaks (i.e.
260 shoulder, stutter or secondary peaks bordering the authentic peak in a same individual) lead
261 to the definition of numerous extra-bins. In such a situation, artefactual peaks can cause the
262 local definition of extra bins into which homologous peaks can be inconsistently assigned.
263 We advise to use MinBin values ranging between 1 and 1.5 bp. The binning is launched
264 through the graphical interface (menu “RawGeno/2. Scoring”) or using the following
265 command lines:

266 Proceed to binning

267 `EXTRACTAFLP(all.dat=AFLP$all.dat,`
268 `samples.names=AFLP$samples.names, MAXBIN=2, MINBIN=1)`

269

270 View results (that are assigned into a “data.binary” object)

271 `attributes(data.binary)`

272 `data.binary$data.binary`

273

274 Binning is an automated and straightforward analysis step that users might want to review

275 interactively. RawGeno includes a visualization device for manually editing the binning by adding,

276 removing or modifying the width and position of bins (this tool is only available from the graphical

277 interface, Fig 2 b). This device includes several help-to-decision statistics such as the average size

278 and the number of presences associated with each bin.

279

280 *2.4.5. Filtering*

281 Once defined, bins can be filtered according to their properties and/or quality. Note that such

282 filtering strategies require analyzing AFLP reactions with a consistent quality across individuals.

283 Filtering options are accessible from the scoring menu when using the graphical user interface

284 (menu “RawGeno/2. Scoring”). Command line users will set accordingly the EXTRACTAFLP

285 function.

286 Proceed to binning and filtering simultaneously

287 `EXTRACTAFLP(all.dat=AFLP$all.dat,`

288 `samples.names=AFLP$samples.names, MAXBIN=2, MINBIN=1, RMIN=100,`

289 `RMAX=500, cutRFU=50, who='B', thresh=95)`

290

291 In its current version, RawGeno includes three kinds of filters.

292 a. The size filter restricts binning to a given portion of the electropherogram (RMIN and

293 RMAX parameters). We advise to limit the binning to peaks included in the range of the size

294 ladder, because their size is accurately interpolated by PeakScanner (in contrast to larger

295 peaks where the size is extrapolated). We recommend discarding peaks with small sizes (i.e.

296 smaller than 100 bp, RMIN=100) as they are more likely to be homoplastic **(9, 10)**. In

addition, large size peaks should as well be considered cautiously because their fluorescence intensity might not always be consistent across individuals. Because this upper limit might vary according to datasets, we advise to run preliminary analyses and check electropherograms to confidently determine it.

- b. The second filter eliminates bins according to their average fluorescence (i.e. the cutRFU parameter). This filter assumes that bins with a high average fluorescence retrieve a more consistent signal than bins with a low fluorescence. The rationale for this strategy is the following. The fluorescence of an AFLP fragment largely determines its detection probability during the PeakScanner analysis of electropherograms. Therefore, fragments that systematically produce low fluorescences are more likely to be erroneously recorded as absent from electropherograms as they might pass the threshold in some reactions but not in others just by chance. The computation starts by normalizing fluorescence intensities across samples using the sum normalization method (6), before computing the average fluorescence of each bin. Hence, keep in mind that the cutRFU parameter applies on normalized values and does not scale with fluorescence measures provided in PeakScanner. Setting this filter is dataset-dependent and we recommend running several trials before producing a definitive dataset (see Note 4). Refer to the works of Whitlock *et al.* (6) and Herrmann *et al.* (7) for more sophisticated filtering R scripts based on fluorescence intensities. Bridging RawGeno with these algorithms is achieved by exporting fluorescence results instead of a binary matrix. Once binning is achieved, use the following command lines.

Retrieve raw fluorescence results, stored in a matrix corresponding to the usual binary matrix

```
mat.rfu=t(data.binary$data.height.raw)
```

For normalised fluorescences (sum normalization), use instead

```
mat.rfu=t(data.binary$data.height)
```

324 Prepare for export and save as a text-tabulated file (refer to programs' documentation for
325 properly formatting files).

```
326 mat.rfu[is.na(mat.rfu)==T]=0  
327 write.table(mat.rfu, "MyFluorescenceFile.txt", quote=F,  
328 sep="\t")
```

329

330 c. The reproducibility filter evaluates bin quality according to their robustness across AFLP
331 reactions, by relying on replicated samples. This filter assumes that replicated individuals
332 were selected randomly from the original dataset, so as to scan the genetic diversity at best
333 (see Notes 2 and 4). Keep in mind that RawGeno identifies replicated individuals using their
334 names. Replicates must be named using the original individual name plus a suffix letter. The
335 suffix is matched using the “who” parameter of the filtering algorithm. As an example,
336 “mysample.fsa” and “mysampleB.fsa” are a pair of original-replicated samples, being
337 identified with a “B” suffix (therefore, set who = “B” when filtering). For each bin, RawGeno
338 compares original to replicated individuals and calculates the percentage of original-
339 replicated pairs for which the AFLP signal is successfully reproduced. Bins where
340 reproducibility cannot reach a satisfactory rate (i.e. the “thresh” parameter, a user-defined
341 reproducibility percentage) are eliminated from the final dataset.

342

343 *2.4.6. Review of results*

344 RawGeno offers two displays for exploring scoring results (menu “RawGeno/3. Quality
345 Check/Samples Checking”). The binary matrix can be directly visualized using a heatmap, showing
346 individuals sorted according to their genetic relatedness. Alternatively, individuals can be examined
347 with a principal coordinates analysis. Both displays allow plotting quality statistics or external
348 information (i.e. picked from the “info table” cited above) against the AFLP results. These displays
349 are only available from the graphical user interface (Fig 3 a, b).

350

2.4.7. Exporting files

RawGeno includes functions for producing binary tables and standard exports for Arlequin, Hickory, Popgen, AFLPsurv, STRUCTURE 2.2, Mltr, Spagedi, Dfdist, Treecon, Baps, PAUP, Structurama, MrBayes and NewHybrids **(11)**. Furthermore, these exports can be sliced according to information provided in the info table and produce ad-hoc subsets. These functions are accessible from the graphical interface (“RawGeno/4. Save”) or using the following command lines.

Retrieve the binary matrix from RawGeno and import the info table

```
mat01=t(data.binary$data.binary)
```

```
matinfo=read.delim("MyInfoTable.txt", header=T)
```

Cross-reference the AFLP results to the info table

```
popsA = row.names(mat01)
```

```
popsB = matinfo$Tag
```

```
mat01= mat01[match(intersect(popsA, popsB), popsA), ]
```

```
matinfo = matinfo[match(intersect(popsA, popsB), popsB), ]
```

Remove monomorphic bins

```
mat01=mat01[ , colSums(mat01)>0 & colSums(mat01)<nrow(mat01)]
```

Produce the required outputs, e.g. for STRUCTURE 2.2. (refer to the library documentation for further details regarding exporting functions).

```
Structure.popsD(mat01, pops=matinfo$MyPopsColumn, path=getwd(),
```

```
name="MyStructure2.2File.txt")
```

Users willing to analyse AFLP signals as codominant markers **(12)** should use command lines described above to export fluorescence data (a proxy of allele copy number in genomes) associated with binary matrices.

379 2.4.8. Handling data from previously scored projects

380 Users willing to merge, visualize and/or produce exports from datasets that were already scored
381 can import binary tables within RawGeno using the “RawGeno/1. Files/Import” menu. From the
382 command line, such an operation is done as follows.

383 Select files to merge

```
384 list.merge=tk_choose.files(caption='Choose Files to Merge')
```

385

386 Or specify a directory in which the binary matrices are stored.

```
387 mypath="C:/MyDocuments/MyBinaryMatricesDirectory"
```

```
388 list.merge=dir(mypath, pattern='.txt')
```

389

390 Proceed to merging

```
391 MERGING(transpose = "indRows", exclude = T, replacewith = NA)
```

392

393 The transpose parameter states whether the binary matrices store individuals as lines
394 (“indRows”) or as columns (“indColumns”), the exclude parameter defines whether individuals that
395 are not shared by all matrices will be removed from the final merged dataset (exclude = “T”). If
396 kept (exclude = “F”), individuals with missing AFLP genotypes will be completed using NA values
397 (replacewith = NA) when no data is available. Note that the merged matrix is stored into a
398 “mergedTable” object. Visualization and exports can be performed using the graphical user
399 interface, as described above. Command lines for exporting merged matrices are given below.

```
400 mat01=mergedTable
```

```
401 matinfo=read.table("MyInfoTable.txt", header=T)
```

```
402 popsA = row.names(mat01)
```

```
403 popsB = matinfo$Tag
```

```
404 mat01= mat01[match(intersect(popsA, popsB), popsA), ]
```

```
405 matinfo = matinfo[match(intersect(popsA, popsB), popsB), ]
```

```
406 mat01=mat01[ , colSums(mat01)&colSums(mat01)<nrow(mat01)]
```



```
407     Structure.popsD(mat01, pops=matinfo$MyPopsColumn, path=getwd(),  
408 name="MyStructure2.2File.txt")
```

409

410 **3. Conclusions and perspectives**

411 We present here a complete suite of tools to automate the scoring of AFLP datasets using free
412 software applications. Our program proposes an integrated solution to manage all the components of
413 the analysis pipeline. Accordingly, samples are checked by removing non-satisfactory
414 electropherograms at the very beginning of the analyses and AFLP genotypes are associated with
415 user-specified information while producing ad-hoc exports. Bins are managed using an automated
416 algorithm and can be edited manually using a dedicated graphical user interface.

417 As a next milestone, we plan to develop RawGeno into two complementary directions:
418 incorporating the handling of electropherograms (which is now is part of PeakScanner) and the
419 enhancement of bin filtering possibilities. Indeed, the RawGeno version currently under
420 development already integrates functions for detecting and calculating the size of AFLP peaks along
421 electropherograms. Finally, achieving connections with the R scripts of Herrmann *et al.* (6) and
422 Whitlock *et al.* (7) is another way to improve RawGeno.

423

4. Notes

Note 1

Linux users might need to run R as “sudo” users to properly install companion packages (i.e. *vegan* and *tkrplot*). In addition, troubles might arise because R libraries are downloaded as source code and compiled locally before being installed. This requires that all compilers needed by R (such as *gc*, *gcc*, *gcc-fortran* and others) have been installed locally, before attempting the installation of external R packages. In OpenSUSE, the necessary compilers can be obtained using YaST2 (into the *rpm* groups dedicated to development tools). Ubuntu users are more fortunate because Synaptic Manager can install ready-to-use R libraries in addition to usual compilers (refer to <http://cran.r-project.org/bin/linux/ubuntu/README> for further details regarding repository addresses).

Note 2

Whereas manual scoring allows permanent but subjective adjusting of the criteria defining whether an AFLP peak should be recorded as present or absent, automated peak detection algorithms apply uniform fluorescence thresholds. This requires datasets showing low AFLP quality variation, because fluorescence differences between samples will be reflected in the final binary matrix. For instance, we encountered problematic situations with AFLP reactions showing different fluorescence offsets among PCR plates. These “plate effects” can be highlighted with principal coordinate analyses, as samples are clustered according to PCR plates (i.e. due to plate-specific losses of AFLP bands). These situations are especially difficult to handle without repeating wet-lab experiments. Solving this problem computationally remains difficult and asks for setting sample (or plate)-specific detection sensitivities, which is beyond PeakScanner possibilities. The present version of RawGeno conservatively proposes the removal of samples with unusual AFLP profiles, by relying on the distribution of the number of AFLP peaks. Future developments of RawGeno will attempt analysing electropherograms directly. The problem highlights the crucial importance of using robust and standardized lab protocols. In the following, we list several tips helping to retrieve consistent results from AFLP reactions, when using automated scoring [also consult Bonin *et al.*

451 (13) and Gugerli *et al.* (14)].

- 452 a. Randomize reactions on PCR plates, in order to properly discriminate technical bias from
453 biological signals.
- 454 b. Standardize all reaction steps: adjust DNA concentrations after spectrometer quantification,
455 limit impacts of pipetting errors by preparing reaction mixes in batches, run critical reactions
456 in uniform conditions (for instance perform restriction steps in an incubator rather than in a
457 thermocycler) and, importantly, run PCRs on the same thermocycler.
- 458 c. Optimize signal detection during genotyping analysis, for instance by increasing the injection
459 time of automated sequencers at the beginning of electrophoresis, and prefer primer pairs
460 showing strong and consistent amplifications.

461

462 **Note 3**

463 Based on empirical case-study datasets (available from the authors upon request), we provide
464 help-to-decision statistics aimed at guiding users in setting their RawGeno analysis. We analysed 17
465 AFLP primer datasets: *Aegilops geniculata*, two primer pairs (5); *Arum* spp., two primer pairs
466 (Espindola *et al.* unpubl. data); *Baldellia* spp., two primer pairs (Arrigo *et al.*, unpubl. data);
467 *Bupleurum ranunculoides*, three primer pairs (Labhardt *et al.*, unpubl. data); *Deschampsia littoralis*
468 and *Deschampsia caespitosa*, each with three primer pairs (15); *Peucedanum ostruthium*, two
469 primer pairs (Borer *et al.*, unpubl. data). These studies included between 87 and 509 individuals
470 (mean: 255) of either intra- and interspecific sampling; 4% to 51% (mean: 23%) of the individuals
471 were replicated. We varied binning parameters (i.e. MinBin and MaxBin) and measured their effects
472 on the width of bins, the datasets' polymorphism (i.e. the proportion of bins with presence
473 frequencies ranging between 5% and 95%) and the bin reproducibility (i.e. the proportion of
474 replicated sample pairs over which the focal bin is successfully reproduced; this measure is
475 independent of the total number of bins in the dataset and is therefore suitable in the context of
476 binning optimization).

477 MinBin and MaxBin acted as boundaries on the width of bins. They determined how accurately

478 an AFLP signal was reflected into the final presence/absence matrix. Both parameters were
479 explored for values ranging between 0.1 bp and 5 bp, with 0.1 bp increments, therefore totalizing
480 1176 “binning trials” per dataset (Fig 4. a, b).

481 Using exaggeratedly small MaxBin values forces the binning algorithm to define narrow bins.
482 This situation leads to “oversplitting”, a bias where AFLP signals are coded into more bins than
483 needed. In this case, an AFLP locus is coded using several adjacent bins appearing as inconsistent
484 when considered independently. Our results showed oversplitting evidence for bin widths below
485 0.5 bp, with decreased bins' polymorphism and reproducibility (Fig 4. c, d). This situation should be
486 avoided, and we recommend using values larger than 0.5 bp for MaxBin (in contrast, the MinBin
487 parameter had little effects on oversplitting).

488 Using exaggeratedly large MinBin and MaxBin values introduces “technical homoplasy” [as
489 defined in (4)], a bias where AFLP signals are coded using less bins than required. Although
490 merging artefactual secondary peaks with authentic peaks is desirable, a process that increases the
491 consistency of binning (see above), exaggerated merging tends to artificially increase similarity
492 between unrelated samples and has immediate effects on AFLP polymorphism and reproducibility.
493 Our results showed that technical homoplasy was reflected by a decrease in bin polymorphism,
494 when MinBin and MaxBin both exceeded 2 bp (i.e. corresponding to bin widths larger than 0.8 bp).
495 On the other hand, reproducibility increased along with technical homoplasy due to the addition of
496 some level of artefactual similarity among samples.

497 From these results, we suggest to screen binning parameters by considering bin polymorphism as
498 a main optimization criterion (Fig 4. a, c). Reproducibility statistics should not be considered for
499 binning optimization because of their inability to detect technical homoplasy (Fig 4. b, d). See
500 Holland *et al.* (7) for further considerations about binning optimization.

501

502 **Note 4**

503 In its current version, RawGeno includes three filters that can be applied after binning has been
504 achieved. These are a size (i.e. the region of electropherograms over which the analysis must be

505 carried out), a fluorescence and a reproducibility filters. We applied filters to the 17 datasets
506 explored during binning optimizations (see Note 3). The size filter was not tested and all datasets
507 were analysed for bins ranging between 100 and 400 bp (i.e. corresponding to the electrophoresis
508 region where size interpolation is accurate). The two remaining filters were explored starting from
509 binned datasets that maximized polymorphism.

510 The first filter considers bins' quality to vary according to average fluorescence. Indeed, bins
511 with an average fluorescence close to the detection threshold used in PeakScanner are more likely
512 to reflect inconsistent signals (e.g. technical false absences) than bins with strong average
513 fluorescence (see above). The filter removes bins according to a user-defined lower fluorescence
514 limit. Keep in mind that this limit applies on normalized values. We tested limits ranging between 0
515 (i.e. no filtering) and 400 rfu (with 10 rfu increments) and measured polymorphism and
516 reproducibility changes caused by this filtering (Fig 5. a, b). All but one dataset gained in
517 polymorphism and reproducibility when filtering was optimized. However, while filtering increased
518 dataset reproducibility, it drastically decreased polymorphism when set up improperly. Its use hence
519 requires dataset-depending optimization to limit information reduction in datasets. We recommend
520 testing fluorescence thresholds below 200 rfu.

521 The reproducibility filter assesses the robustness of AFLP signals and removes bins that are not
522 satisfactorily reproducible. We tested filtering using bin reproducibility limits ranging between 0
523 and 100% (with 5% increments) and measured polymorphism and global error rate of datasets
524 (Fig 5. c, d). We present error rates (**13**) instead of bin reproducibility here because error rates are
525 more commonly reported in AFLP studies. Filtering effectively reduced the global error rate along
526 with the removal of non-reproducible bins from datasets. As expected, filtering also decreased the
527 polymorphism of datasets. Nevertheless, this filter generally removed numerous bins, and we
528 clearly propose avoiding the use of reproducibility values larger than 95%. Our experience shows
529 that using 85% as a threshold provides satisfactory results by balancing the error rates of datasets
530 (i.e. 5% in average, which we consider reasonable) with polymorphism. A more stringent filtering
531 can be applied (e.g. in studies where peak reproducibility is especially relevant such as in genome

532 scans), but it might affect the information content of datasets and thus requires the inclusion of
533 additional primer pairs to achieve a reasonable level of polymorphism.

534

535 **5. Acknowledgments**

536 We are grateful to Anahi Espindola, Matthias Borer and Aurélien Labhardt for kindly providing
537 the AFLP datasets presented in Notes 2-4. We thank Nicolas Salamin for having provided a working
538 environment to N. Arrigo during the redaction of this manuscript. Felix Gugerli, Sabine Brodbeck,
539 Piya Kittipadakul, Julie B. Hébert, Julie Lee Yaw, Christian Parisod, Roland Dubillard and Robin
540 Arnoux provided helpful comments during the development of RawGeno. This project was
541 supported by the National Centre of Competence in Research (NCCR) Plant Survival and the PNR-
542 59, two programs of the Swiss National Science Foundation.

6. References

1. Vos P, Hogers R, Bleeker M, Reijans M, van de Lee T, Hornes M, Frijters A, Pot J, Peleman J, Kuiper M, Zabeau M (1995) AFLP: a new technique for DNA fingerprinting. *Nucleic Acids Res* 23:4407–4414.
2. Pompanon F, Bonin A, Bellemain E, Taberlet P (2005) Genotyping errors: causes, consequences and solutions. *Nat Rev Genet* 6:847–859.
3. Meudt HM, Clarke AC (2007) Almost forgotten or latest practice? AFLP applications, analyses and advances. *Trends Plant Sci* 12(3):106-117.
4. Arrigo N, Tuszynski JW, Ehrich D, Gerdes T, Alvarez N (2009) Evaluating the impact of scoring parameters on the structure of intra-specific genetic variation using RawGeno, an R package for automating AFLP scoring. *BMC Bioinformatics* doi:10.1186/1471-2105-10-33.
5. Arrigo N, Felber F, Parisod C, Buerki S, Alvarez N, David J, Guadagnuolo R (in press) Origin and expansion of the allotetraploid *Aegilops geniculata*, a wild relative of wheat. *New Phytol*.
6. Whitlock R, Hipperson H, Mannarelli M, Butlin RK, Burke T (2008) An objective, rapid and reproducible method for scoring AFLP peak-height data that minimizes genotyping error. *Mol Ecol Res* 8:725-735.
7. Herrmann D, Poncet BN, Manel S, Rioux D, Gielly L, Taberlet P, Gugerli F (2010) Selection criteria for scoring amplified fragment length polymorphisms (AFLPs) positively affect the reliability of population genetic parameter estimates. *Genome* 53(4):302-310.
8. Holland BR, Clarke AC, Meudt HM (2008) Optimizing automated AFLP scoring parameters to improve phylogenetic resolution. *Syst Biol* 57(3):347-366.
9. Vekemans X, Beauwens T, Lemaire M, Roldan-Ruiz I (2002) Data from amplified fragment length polymorphism (AFLP) markers show indication of size homoplasy and of a relationship between degree of homoplasy and fragment size. *Mol Ecol* 11(1):139-151.
10. Paris M, Bonnes B, Ficetola GF, Poncet BN, Després L (2010) Amplified fragment length

569 homoplasmy: in silico analysis for model and non-model species. BMC Genomics
570 doi:10.1186/1471-2164-11-287.

571 11. Ehrlich D (2006) AFLPdat: a collection of R functions for convenient handling of AFLP
572 data. Mol Ecol Notes 6: 603-604.

573 12. Gort G, van Eeuwijk F (2010) Codominant scoring of AFLP in association panels. Theor
574 Appl Genet doi:10.1007/s00122-010-1313-x.

575 13. Bonin A, Bellemain E, Eidesen PB, Pompanon F, Brochmann C, Taberlet P (2004) How to
576 track and assess genotyping errors in population genetics studies. Mol Ecol 13:3261–3273.

577 14. Gugerli F, Englisch T, Niklfeld H, Tribsch A, Mirek Z, Ronikier M, Zimmermann NE,
578 Holderegger R, Taberlet P, IntraBioDiv Consortium (2008) Relationships among levels of
579 biodiversity and the relevance of intraspecific diversity in conservation - a project synopsis.
580 Perspect Plant Ecol Evol Syst 10: 259-281.

581 15. Peintiger M, Arrigo N, Brodbeck S, Koller A, Imsand M Holderegger R (2010) Genetische
582 und morphologische Differenzierung der endemischen Grasart *Deschampsia littoralis*
583 (Gaudin) Reut. – Wie verschieden sind die Population am Bodensee und am Lac de Joux im
584 Vergleich zu *D. cespitosa* (L.) P. Beauv.? Vogelwarte Radolfzell, University of Neuchâtel,
585 WSL.

7. Figure captions

Fig 1. Checking quality of AFLP reactions. (a) Interactive device for filtering samples, before performing binning. Samples are sorted according to the number of AFLP peaks successfully amplified. Summary statistics (i.e. 5%, 50% and 95% quantiles) are provided to help users selecting samples to be included in further analyses. (b) Visualization of 96-well PCR plates. Four quality statistics are available in RawGeno (number of AFLP peaks per sample, mean and variance of fluorescence intensity and outlier detection index; see RawGeno manual). Boxplots provide summary statistics for rows and columns of PCR plates, respectively.


Fig 2. Binning algorithm implemented in RawGeno (from (4)). (a) Left panel: main steps followed by the algorithm to define bins; right panel: illustration of binning with two samples (S1 and S2); the bin widths (i.e. the difference in size between the longest and the shortest amplicons included in the considered bin) and the technical homoplasmy rates (i.e. HR, the mean number of peaks belonging to the same sample that are included in a same bin) are indicated. (b) Binning edition device, allowing users to translate, resize, add or remove bins interactively, starting from bins that were initially defined by the automated algorithm. Summary statistics (i.e. the average size of bins and the number of AFLP peaks present per bin) are provided as editing guidelines.

Fig 3. Reviewing results. RawGeno includes basic visualization devices for performing preliminary data mining. Specifically, results can be reviewed using (a) heatmaps of the binary matrix, where samples are sorted according to their genetic similarity and (b) principal coordinates analysis of the corresponding matrix. Both devices can compare AFLP results with either quality statistics (i.e. number of AFLP peaks per sample, mean and variance in fluorescence intensity and outlier detection index) or external information provided by users (e.g. the population from where samples were collected). In addition, both devices are handled through a graphical user interface for sorting and selecting samples to be visualized.


Fig 4. Binning parameters. We analysed 17 different datasets (see Note 3) in RawGeno, by varying binning parameters that have an effect on the width of bins (i.e. MinBin and MaxBin) and measured the associated effects on bin width, polymorphism and reproducibility. All results were corrected in order to normalize statistics among datasets (i.e. using non-linear mixed effect models, with the dataset origin considered as a covariable). Upper row: triangle plots displaying (a) bin polymorphism and (b) reproducibility averaged for the 17 datasets as a function of MinBin and MaxBin parameters. Lower row: scatterplots of corresponding (c) bin polymorphism and (d) bin reproducibility statistics, displayed according to the median bin width of binning trials. 5% and 95% confidence intervals (dashed lines) and the average (continuous line) are displayed. (a) to (d) Trials optimizing bin polymorphism are represented as dots for each of the 17 datasets (boxplots indicate bin widths associated to optimized trials).

Fig 5. Filtering parameters. We filtered datasets (using optimal binning parameters as defined in Note 3) by increasing bin fluorescence and reproducibility thresholds. Median bin polymorphism and error rates were measured for the final datasets (dataset-specific offsets were not corrected). Upper row: bins are filtered according to their average fluorescence. Results are displayed as a comparison with non-filtered trials, with “D bin polymorphism” and “D bin reproducibility” being the difference between filtered and non-filtered datasets. Therefore, values greater than zero reflect cases where filtering increased the median (a) polymorphism or (b) reproducibility of bins in the final dataset. Conversely, values lower than zero indicate decreases in these statistics associated with filtering. Lower row: filtered bins according to their reproducibility and consequences on median bin (c) polymorphism and (d) error rate of the final datasets. Results of each dataset (grey lines) are displayed along with summary statistics (dashed lines: 5% and 95% confidence intervals; continuous line: average). (a) to (d) Trials optimizing bin polymorphism after filtering are presented as dots for each of the 17 datasets (boxplots indicate fluorescence threshold associated to optimized filtering).


(a)




(b)




(a)




(b)




(a) Median bin polymorphism



(b) Median bin reproducibility



(c)



(d)

